STANDARD Q**COVID-19 IgM/IgG Duo**

STANDARD™ Q COVID-19 IgM/IgG Duo Test

PLEASE READ CAREFULLY BEFORE YOU PERFORM THE TEST

SD BIOSENSOR


REF Q-NCOV-01D

KIT CONTENTS**PREPARATION - Be sure to test both STANDARD Q COVID-19 IgM and IgG simultaneously.**

- 1** Carefully read instructions for using the STANDARD Q COVID-19 IgM/IgG Duo Test.

- 2** Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.


- 3** Open both STANDARD Q COVID-19 IgM and IgG pouches, and check the test devices and the desiccant in each pouches.

**TEST PROCEDURE - Be sure to test both STANDARD Q COVID-19 IgM and IgG simultaneously.**


The test procedures for both COVID-19 IgM and IgG are the same.

Using Capillary whole blood


- 1** **Collecting of Specimen**
Using a capillary tube, collect the 10μl of capillary whole blood to the black line of the capillary tube.

- 2** **Adding of Specimen**
Add the collected capillary whole blood to the specimen well of the test device.


- 3** **Dropping of buffer**
Add 3 drops (90μl) of buffer into the buffer well of the test device.


- 4** **Reading Time**
Read test result at 10~15 minutes.

**Using serum/plasma/venous whole blood**


- 1** **Collecting of Specimen**
Using a micropipette, collect the 10μl of serum, plasma or venous whole blood.



- 2** **Adding of Specimen**
Add the collected serum/plasma or venous whole blood to the specimen well of the test device.


- 3** **Dropping of buffer**
Add 3 drops (90μl) of buffer into the buffer well of the test device.

- 4** **Reading Time**
Read test result at 10~15 minutes.


**INTERPRETATION OF TEST RESULT****Positive**

IgM Positive

Negative

IgG Positive

IgM Negative

Invalid

IgG Negative

IgM Positive

Re-test with a new test device.

1. A colored band will appear in the top section of the result window to show that the test is working properly. This band is control line (C).
2. A colored band will appear in the lower section of the result window. These bands are test line of IgM/IgG (M, G).

3. Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.

- * STANDARD Q COVID-19 IgM/IgG Duo Test may cross-react with antibody against SARS-CoV-1.

- * Results from antibody testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform infection status.

- * Positive results should be considered in conjunction with the clinical history, RT-PCR results and other data available.

EXPLANATION AND SUMMARY**[Introduction]**

Coronavirus is a single-stranded positive-sense RNA virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen in humans, animals, pets, and human disease. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or "COVID-19", was discovered due to Wuhan Viral Pneumonia cases in 2019 and was named by the World Health Organization on January 12, 2020. WHO confirmed that COVID-19 can cause colds, the Middle East Respiratory Syndrome (MERS) and more serious diseases such as severe acute respiratory syndrome (SARS). This kit is helpful for the auxiliary diagnosis of coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

[Intended use]

STANDARD Q COVID-19 IgM/IgG Duo Test is a rapid chromatographic immunoassay for the qualitative detection of specific antibodies to SARS-CoV-2 present in human serum, plasma or whole blood. This test is for in vitro professional diagnostic use and intended as an aid to diagnosis of SARS-CoV-2 infection in convalescent phase of patient with clinical symptoms with SARS-CoV-2 infection. It provides only an initial screening test. More specific alternative diagnosis methods should be performed in order to obtain the test line.

[Test principle]

STANDARD Q COVID-19 IgM/IgG Duo Test has two pre-coated lines, "C" Control line, "G" Test line for the COVID-19 IgG Device and "C" Control line, "M" Test line for the COVID-19 IgM device on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimens. Monoclonal anti-COVID-19 antibody is coated on the control line region and Monoclonal anti-human IgG is coated on the test line region in IgG device and Monoclonal anti-human IgM is coated on the test line region in IgM device. Recombinant COVID-19 nucleocapsid protein conjugated with colloidal gold particles are used as detectors for COVID-19 IgG Device and COVID-19 IgM Device. During the test, SARS-CoV-2 antibodies in the specimen react with recombinant COVID-19 nucleocapsid protein with colloidal gold particles making antibody-antigen gold particle complex. This complex migrates on the membrane via capillary action along the test line, where it will be captured by the Monoclonal anti-human IgG antibody or Monoclonal anti-human IgM antibody. A violet test line would be visible in the result window if SARS-CoV-2 antibodies are present in the specimen. The intensity of violet test line will vary depending upon the amount of SARS-CoV-2 antibodies present in the specimen. If SARS-CoV-2 antibodies are not present in the specimen, no color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

[Kit contents]

① Test device (individually in a foil pouch with desiccant) ② Buffer bottle ③ Capillary tube (10μl) ④ Instructions for use

KIT STORAGE AND STABILITY

Store the kit at room temperature, 2-30°C / 36-86°F, out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

WARNINGS AND PRECAUTIONS

1. Do not reuse the test kit.
2. Do not use the test kit if the pouch is damaged or the seal is broken.
3. Do not use the buffer of another lot.
4. Do not smoke, drink or eat while handling specimen.
5. Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Handle the test kit carefully to avoid damage to the test device.
8. Observe established precautions against microbiological hazards throughout testing procedures.
9. Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
10. Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
11. Good laboratory practice recommends the use of the control materials. Users should follow the appropriate federal state, and local guidelines concerning the frequency of assaying external quality control materials.

SPECIMEN COLLECTION AND PREPARATION**[Serum]**

1. Collect the whole blood into the commercially available plain tube. NOT containing anti-coagulants such as heparin, EDTA, Sodium citrate by venipuncture and leave to settle for 30 minutes for clot retraction and then centrifuge blood to get serum specimen of supernatant.
2. If serum in the plain tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 1 week after collection. Using the specimen in the long-term keeping more than 1 week can cause non-specific reaction. For prolonged storage, it should be at below -40°C/-40°F.
3. They should be brought to room temperature prior to use.

[Plasma]

1. Collect the venous blood into the commercially available anti-coagulant tube such as heparin, EDTA, Sodium citrate by venipuncture.
2. If plasma in an anti-coagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 1 week after collection. Using the specimen in the long-term keeping more than 1 week can cause non-specific reaction. For prolonged storage, it should be at below -40°C/-40°F.
3. They should be brought to room temperature prior to use.

[Whole blood]

- **Capillary whole blood**
 1. Capillary whole blood should be collected aseptically by fingertip.
 2. Clean the area around the fingertip with alcohol swab.
 3. Squeeze the end of the fingertip and pierce with a sterile lancet.
 4. Using a capillary tube, collect the 10μl of capillary whole blood to the black line of the capillary tube.
 5. The capillary whole blood must be tested immediately after collection.
- **Venous whole blood**
 1. Collect the venous blood into the commercially available anti-coagulant tube such as heparin, EDTA, Sodium citrate by venipuncture.
 2. If venous whole blood in an anti-coagulant tube is stored in a refrigerator at 2-8°C/36-46°F, the specimen can be used for testing within 1-2 days after collection.
 3. Do not use hemolyzed blood Specimens.

- As known relevant interference, hemolytic Specimen, rheumatoid factors-contained Specimen and lipemic Specimen can lead to impair the test results.
- Separate disposable materials for each Specimen in order to avoid cross-contamination which can cause erroneous results.

PERFORMANCE CHARACTERISTICS**[Clinical evaluation]**

Performance characteristic for the STANDARD Q COVID-19 IgM/IgG Duo Test for rapid detection of anti-SARS-CoV-2 antibodies was established in retrospective, multi-center, randomized, single-blinded study conducted at a trial site in KOREA during the 2020 SARS-CoV-2 pandemic situation. A total of 504 retrospective specimens were tested using the STANDARD Q COVID-19 IgM/IgG Duo Test. These specimens consisted of serum from PCR positive or negative confirmed patients. The performance of the STANDARD Q COVID-19 IgM/IgG Duo Test was compared to a commercialized molecular assay. Although the STANDARD Q COVID-19 IgM/IgG Duo Test allows to test for IgM and IgG separately, due to the differing inter-patient time response to the virus, any individual with positive result for the IgM or the IgG test should be read as a positive for anti-SARS-CoV-2 antibodies. The combined test result (positive for IgM and/or IgG or negative for IgM and/or IgG) was used to calculate the total test sensitivity and specificity.

[Test sensitivity]

The seroconversion time of IgM and IgG antibodies varies from person to person, but it was estimated to be around 7 days after onset of symptom¹. The STANDARD Q COVID-19 IgM/IgG Duo Test showed 94.33% of sensitivity using specimens from patients 7 days after symptom onset (combined IgM+IgG).

Table 1. Summary of the sensitivity of the STANDARD Q COVID-19 IgM/IgG Duo Test compared to PCR using specimens 7 days after symptom onset.

	PCR		
	Positive	Negative	Total
STANDARD Q COVID-19 IgM/IgG Duo Test	183	0	183
Positive	11	0	11
Total	194	0	194
94.33%			(183/194, 95% CI, 90.08% - 97.14%)

Table 2. Summary of the sensitivity of the STANDARD Q COVID-19 IgM/IgG Duo Test compared to PCR using specimens with less than 7 days after symptom onset.

	PCR		
	Positive	Negative	Total
STANDARD Q COVID-19 IgM/IgG Duo Test	31	0	31
Positive	14	0	14
Total	45	0	45
68.89%			(31/45, 95% CI, 53.35% - 81.83%)

Table 3. Summary of the sensitivity of the STANDARD Q COVID-19 IgM/IgG Duo Test compared to PCR using specimens in the period of 7 days to 14 days after symptom onset.

	PCR		
	Positive	Negative	Total
STANDARD Q COVID-19 IgM/IgG Duo Test	73	0	73
Positive	10	0	10
Total	83	0	83
87.95%			(73/83, 95% CI, 78.96% - 94.07%)

Table 4. Summary of the sensitivity of the STANDARD Q COVID-19 IgM/IgG Duo Test compared to PCR using specimens 14 days after symptom onset.

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